

# **Exhibit A**


# RECOVERY<sup>†</sup> G2 FILTER SYSTEM

Timeless Performance<sup>†</sup>

## NOW AVAILABLE - JUGULAR DELIVERY SYSTEM



The Recovery<sup>†</sup> G2 Filter combines **the best design features** of Bard's existing vena cava filters to create a **brand-new filter platform** — taking strength and stability to a new level.



The newly enhanced Recovery<sup>†</sup> G2 Filter continues the Bard tradition of filter innovation spanning over a decade.

- Increased **MIGRATION RESISTANCE\***
- **IMPROVED CENTERING\***
- Enhanced **FRACTURE RESISTANCE\***

\* Data on File

### Clot Trapping and Caval Patency

The Recovery<sup>†</sup> G2 Filter utilizes the proven conical filter shape arranged into two offset layers that effectively trap large and small emboli without compromising caval patency.

### Secure Fixation

Now featuring a wider leg span and thicker fixation hooks, the newly enhanced Recovery<sup>†</sup> G2 Filter resists migration across an even broader range of caval distension and higher pressures.\*

\* Maximum indicated caval diameter is 28 mm. Data on File



### Low-Profile

7F delivery system is the lowest profile of any conical filter on the market.

### Self-Centering

Specially designed pusher wire and articulated arms promote a centered filter placement, even through tortuous anatomy.

<sup>†</sup> Bard, Recovery and Timeless Performance are registered trademarks or trademarks of C. R. Bard, Inc. or an affiliate. Copyright 2005 C. R. Bard, Inc. All rights reserved. **Legal Notice**

Please consult product labels and package inserts for indications, contraindications, hazards, warnings, cautions, and information for use.

Copyright © 2006, C. R. Bard, Inc. All Rights Reserved